

MX8 Cardiovascular Ultrasound System

Transducer Family



CW2 (Non-imaging CW)

Applications: Adult Cardiac,
Pediatric Cardiac
CW frequency: 2.0 MHz
Biopsy guide: NA



P10-4s (Phased)

Applications: Adult Cardiac,
Pediatric Cardiac
FOV: 90°
Biopsy guide: NA



P7-3Ts (TEE)

Applications: Adult Cardiac
FOV: 90°
Biopsy guide: NA



SP5-1Ns (Phased)

Applications: Cardiac, Vascular
FOV (max): 90°
Biopsy guide: NGB-011, stainless
steel, undetachable



L9-3s (Linear)

Applications: Vascular
FOV: 43 mm
Convex radius: 20 mm
Biopsy guide: NGB-034, multi-angle,
reusable



L15-3RCs (Linear)

Applications: Vascular
FOV: 38 mm
Biopsy guide: NGB-043, stainless
steel, detachable; NGB-044, stainless
steel, detachable, out-of-phase



P8-2s (Phased)

Applications: Cardiac
FOV: 90°
Biopsy guide: NA



L12-3RCs (Linear)

Applications: Vascular
FOV: 38 mm
Biopsy guide: NGB-043, stainless
steel, detachable; NGB-044, stainless
steel, detachable, out-of-phase



SC5-1Ns (Convex)

Applications: Vascular

FOV: 73 mm

Convex radius: 60 mm

Biopsy guide: NGB-022, stainless steel, detachable

Mindray North America

800 MacArthur Blvd., Mahwah, NJ 07430

Tel: 800.288.2121 Support: 877.913.9663 Fax: 800.926.4275 www.mindray.com

Mindray® is a registered trademark of Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All brands and product names are trademarks of their respective owners. ©2025 Mindray DS USA, Inc. Subject to change. 10/25 P/N: 0002-08-40625 Rev 2.0

mindray
healthcare within reach